Hydro-Gen Fontaine PEM & Inhaler

GEN 3

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GEN 3

Product Manual

Please read this instruction manual carefully before operating this product.

Please keep your product in accordance with the product manual statement.

MIRACLE HAPPENS

For More Information

Scan Here For Details

Follow us :

Visit us :

www.lwgpedia.com | www.livewell-global.com

Safety Precautions

The following precautions are intended to prevent harm or damage to the user and others as well as to the correct use of the product.

- When carrying the Hydro-Gen Fontaine PEM & Inhaler Gen 3 please check the lid and base is tightly closed.
- The jack part of the hydrogen generator (base) is not waterproof.
-Please do not submerge into water.
- Check for any water leakage during operation of the product.
- Do not place heavy objects on top of the product.
- Keep the product out of the reach of children and pet.

WARNING

- Do not use, store or leave in places near fire or where you are exposed to direct sunlight, or where water can be frozen easily.
For example: Bottle sterilizer, inside a hot car, any other place with high temperature or in a refrigerator.
- Do not touch the power plug with wet hands, it may lead to electrical shock.
- Do not use the product if the plug or the wire is damaged.
- Avoid placing the product near magnetic materials.
- Do not disassemble, repair or remodel the product.
- Do not operate the product when there is no water in the container.
- Do not put into dishwasher and bottle sterilizer.
- Do not use metal brushes or scratched sponges when cleaning.

Table of Contents

Benefits of Hydrogen Water and Inhalation Based on Scientific Researches	4
International Certifications of Hydro-Gen Fontaine PEM & Inhaler Gen 3	5
Overview	6
Package Components	6
Directions for Use	7
How to Charge	8
Display	8
How to Use	9
How to Use Inhaler	10
- Benefits of Hydrogen Inhalation	10
How to Use With Bottled Water	10
Settings	
- Volume Setting	11
- Language Setting	11
How to Clean	
- Cleaning the Container	12
- Cleaning of Hydrogen Generator (Base)	12
How to Store	12
Frequently Asked Questions	13


Benefits of Hydrogen Water and Inhalation Based on Scientific Researches

The utilization of hydrogen-rich water is supported by hundreds of leading medical journals. The Japanese and South Korean Health Ministries have approved the use of Hydrogen Water since 1966 and 1978 respectively and this proves that Hydrogen Water has been safe to use for more than 50 years. Molecular Hydrogen Institute (MHI) revealed that hydrogen is proven to be beneficial in more than 200 diseases in almost all human organs in the body. This is supported by 3000 articles and scientific studies involving more than 1600 international researchers.

Prof. Masaru Suzuki, from the Department of Emergency and Critical care Medicine at Keio University Hospital Japan stated that inhalation of Hydrogen Gas can diffuse easily from the lungs to the blood vessels. It is also penetrates the blood brain barrier more easily, so it can be pumped directly throughout the body.

In the 2014 Medical Gas Research journal, a study "A Review of Experimental Studies of Hydrogen As A New Therapeutic Agent in Emergency And Critical Care Medicine", cited recent evidence that methods of administering hydrogen molecules, either through inhalation of hydrogen gas or through drinking hydrogen water, has medical potential as a powerful antioxidant, anti-inflammatory and anti-apoptosis in cells, tissues and organs. Therefore, hydrogen can be used as a therapeutic agent.

International Certifications of Hydro-Gen Fontaine PEM & Inhaler Gen 3

1. FDA (Food and Drug Administration)

A federal agency of the United States Department of Health and Human Service. The FDA is responsible of food safety, supplement, drugs, medical devices, electromagnetic radiation emitting device (ERED), cosmetics, etc.

2. NSF (National Sanitation Foundation)

Certification that states the product have been tested and meet strict standards for global human health protection.

3. KTR (Korea Testing & Research Institute)

Testing and certification institute in South Korea which performs integrate testing, certification, and technical consulting for all fields of the industry.

4. SKYGREEN

Certification for environmentally friendly transparent plastic from SK Chemicals South Korea, a company that develops bio-copolyester plastics that are not harmful to the human body and also reduces global warming.

5. BPA FREE

Certification that states BPA-free products (Bisphenol A), is an industrial chemical that has been used to make certain plastics, such us water bottles or food and beverage containers.

6. LG CHEM

Certification from the largest chemical company in South Korea that provides raw materials, IT and electronic materials and solutions for energy.

7. JOHNSON MATTHEY

Platinum Electrode plate certificate from a multinational chemical and precious metal industry company based in London, England.

8. SHINSHO CORPORATION


Titanium Electrode plates from Japanese company engaged in the export - import and trade of iron, steel, raw materials, nonferrous metals, machinery and welding.

Overview

Package Components

Adapter

USB
Type C
Cable

Bottle
Connector

Nasal Cannula
Connector
(Inhaler Hose Connector)

Directions for Use

The Hydro-Gen Fontaine PEM & Inhaler Gen 3 is packaged with a small amount of distilled water to prevent the hydrogen generator's mesh from drying up and broken. Please do this step if the machine will not be used for more than one day.

WARNING

Please remove the circular plastic cover on top of the hydrogen generating mesh before use. Please keep the circular plastic cover carefully and do not lose it. Only one cover is available for each Hydrogen Fontaine PEM & Inhaler Gen 3 package. The cover is not for sale. Warranty cannot be claimed if the membrane is broken due to drying up.

1. Charging time is between 2-4 hours.
2. Please charge it, when you first use it. (refer 'How to Charge' - page 8.)
3. Please go through one cleaning cycle before first use. (Refer to page 12.)
4. If the product does not work, put a small amount of warm water (40~50°C) covering the mesh part in the container and operate in Mode 1 (3 minutes) 1-2 times or in Mode 2 (10 minutes). Restart the product after generating once.
5. While some water may leak from the oxygen outlet during hydrogen generation, there is no problem. Since it does not mean a malfunction, please dry it and continue use.


WARNING

When operating the product for more than 10 minutes, the pressure inside the container rises, so loosen the container lid to avoid pressure build up. Please don't work continuous. When you open the cap, the water spouts at once and cause machine failure, due to Internal pressure increase.

How to Charge


After a certain period of time, the internal battery will normally be discharged.

Recharging time takes 2-4 hours.

LED battery status will light up while the battery is being charged. When the battery is fully charged, icon will start blinking.

1. To charge the battery, open the rubber jack cover and rotate it upwards to avoid heat during charging (as it is made from environmentally friendly natural rubber)
2. Connect the charger (Adapter and USB cable) into the jack. The battery will light up during charging, then flash when the battery is fully charged.
3. When the battery is fully charged, unplug the USB cable and close the adapter jack.



Adapter Jack USB (USB Type C)


WARNING

Never put metal or liquid inside the USB adapter jack.

Display

Stand by Mode (LED OFF)



Press action button Display will turn on with melody.

How to Use

Mode 1. Generation for 3 minutes

When you want to drink hydrogen water immediately.

1. Press the action button.
The display is turned on and the sound is heard.
2. Press the action button one more time.
The hydrogen generator Mode 1 is activated.



* When you want to stop hydrogen generation,
Press the action button at any time during generation.

Mode 2. Generation for 10 minutes

When you want to inhalation.

1. Press the action button.
The display is turned on and the sound is heard.
2. Press the action button twice more.
The hydrogen generator Mode 2 is activated.


* When you want to stop hydrogen generation,
Press the action button at any time during generation.

How to Use Inhaler

1. Pour mineral water up to the maximum capacity of a bottle.
2. Attach inhaler hose connector on the bottle opening and then put inhaler hose to the nostril.
3. Press the action button.
4. Press the action button twice more.
 - The hydrogen generator Mode 2 is activated.
5. Press again the hydrogen generator Mode 2 for generating until 20 minutes (two times Mode 2).
6. During the electrolysis process, hydrogen gas can be directly inhaled.


Benefits of Hydrogen Inhalation

1. Helps prevent stroke.
2. Helps improve brain function.
3. Helps maintain lung health.
4. Helps reduce digestive disorders.
5. And many others.

How to Use With Bottled Water


1. Connect the bottle connector to the hydrogen generator (base) and to the bottled water.
2. Over turn the bottle such that the hydrogen generator (base) is at the bottom. Press the action button to turn on the power.
3. Press the action button again to activate Mode 1.

Settings

Volume Setting

1. Press the action button to activate it.
2. Press and hold the action button for 3 seconds until the volume icon blinks.
3. Adjust the volume level by pressing the action button to cycle through options.
*Repeat cycle of Mute -> Volume 1 -> Volume 2 -> Volume 3
4. To select a volume, release the action button. It will be automatically set.


Mute -> 1 -> 2 -> 3

Language Setting

1. Press the action button to activate it.
2. Press and hold the action button for 10 seconds until you hear "generating completed" and the battery icon flashes.
3. Set the language by pressing the action button to cycle through options.
*Repeat cycle of English -> Chinese -> Japanese -> Korean.
4. When you select a language, release the action button. It will be automatically set.

English -> Chinese -> Japanese -> Korean

How to Clean

Maintenance is necessary to ensure good performance of the product for a long time. The hydrogen water generator (base) is delicate, clean it carefully.

If minerals start to accumulate in the hydrogen generator (base), it will hinder hydrogen generation and foreign matter and offensive odor may be discharged from the hydrogen generator.

Clean the product according to the following guidelines before starting to use the product.

Cleaning the container

1. Remove the container from the hydrogen generator (base) and wash it with running tap water.
2. If necessary, include neutral detergent in cloth or sponge. Do not use disinfectant or abrasive material.

Cleaning of Hydrogen Generator (base)

1. Fill the bottle with water and add 2-3 drops of citric acid
2. Generate hydrogen in Mode 1 for 1-2 times, then rinse with clean water.
3. Recommended to clean once every 2 weeks, if used daily.

WARNING

Don't soak the product in the water.

How to Store

If you do not use the product for a long time, keep the product upright and put a small amount of water in the product so that the mesh of the hydrogen generator stays wet and close it with the rubber that have been provided.

WARNING

- If the mesh of the hydrogen generator dries out, the product will NOT operate properly.
- Warranty cannot be claimed if the membrane is broken due to drying up.

Frequently Asked Questions

Q. The product does not operate

- A.** Check the battery status of the product.
Please try again after charging the battery. (refer 'How to Charge' - page 8.)

Q. There is water leaking from the bottle

- A.** Check if the hydrogen generator (base) and the body are properly connected.
A small amount of water leak through the oxygen outlet with in the hydrogen generator bottle is acceptable when our product is being used properly. So please continue to use it after drying the bottle.

Q. There is a strange taste in the water, or it smells bad

- A.** Clean the inside of the bottle (refer 'How to Clean' - page 12.)

Q. Hydrogen gas does not come out during production

- A.** The hydrogen generator mesh is dry.
Pour a small amount of lukewarm water (40~50 °C) enough to cover the mesh part in the container and operate 1-2 times in Mode 1. (refer 'Direction for Use' - page 7 point 4.)

Q. There is a slight noise from the product

- A.** Slight noise can flow during the product movement his is normal, not a fault.

It hydrogen generation does not operate efficiently, there is a possibility that minerals may have accumulated inside the generator. Please refer to 'How to Clean' - page 12 and refer to the cleaning method of the hydrogen generator and clean it, before trying again.

Garansi/Warranty*

*Garansi berlaku selama 18 bulan setelah diterima oleh stockist /
Warranty will be validated for 18 months after handed to stockists

Product Specifications

- Hydro-Gen Fontaine PEM & Inhaler
- Power : DC 3.7 V / 1600 mAh (lithium polymer battery)
- Power supply : AC 100-240 V / 50-60 Hz / DC 5 V, 0.3 A
- Dimension : 272 (H) x 75 (D1) x 65 (D2) mm
- Weight : 400 g
- Capacity : 500 ml
- Use temperature : 5-40 °C
- Electrode: Platinum plated + titanium
- 3 minutes/10 minutes : up to 1500 ppb
- Dissolved hydrogen concentration varies with water quality and water temperature.
- KC : MSIP-REM-hca-FT3

For more advanced technical support,

+62 (21) 6570-1843

Head Office of LiveWell Global

Grand Palace, Ruko Palazzo Blok A12
Jl. Benyamin Suaeb, Kav. A5
Kebon Kosong, Kemayoran
Jakarta Pusat, 10630, Indonesia